OPERATING MANUAL
MODE D'EMPLOI
INSTRUCCIONES DE OPERACIÓN

T 814246 04/2014

TruLine ORP 51

TruLine

a xylem brand

TruLine

Operating Manual	3
Technical data	5
Commissioning, measuring, testing	6
Storage	8
Aging	8
Activating platinum electrodes	9
Wear parts and accessories	10
Mode d'emploi	13
Caractéristiques techniques	15
Mise en service, mesure, vérification	16
Stockage	18
Vieillessement	18
Activation d'électrodes de platine	19
Pièces d'usure et accessoires	20
Instrucciones de operación	23
Especificaciones técnicas	25
Puesta en funcionamiento, medición, verificación	26
Almacenamiento	28
Envejecimiento	28
Activar electrodos de platino	29
Piezas de desgaste y accesorios	30
Contact Information	31

TruLine ORP 51

TruLine

ORP ELECTRODE

For the most recent version of the manual, please visit www.yisi.com.

Contact

YSI
1725 Brannum Lane
Yellow Springs, OH 45387 USA
Tel: +1 937-767-7241
800-765-4974
Email: environmental@ysi.com
Internet: www.yisi.com

Copyright

© 2014 Xylem Inc.

Technical data

General data

Model	Reference electrolyte	Electrode material and shape	Junction
TruLine ORP 51	Gel	Platinum pin 1 mm	Fiber

Measurement and application characteristics

Model	Allowed temperature range	Typical application
TruLine ORP 51	-5 ... 80 °C	Field

Shaft dimensions, shaft material, electrical connection

Model	Shaft			Electrical connection		
	Length [mm]	Ø [mm]	Material	Combination electrode connection	Cabel	Cable length
TruLine ORP 51	120	12	Plastic	Plug head	LB1BNC	1 m

Commissioning, measuring, testing

Commissioning

Note

For ORP electrodes with platinum electrode, please observe the instructions in section ACTIVATING PLATINUM ELECTRODES.

Prepare the combination electrode for measuring as follows:

- Remove the watering cap from the electrode tip. Possible salt deposits in the area of the watering cap do not affect the measuring characteristics and can easily be removed with deionized water.

Note

Please keep the watering cap. It is required for storing the combination electrode. Always keep the watering cap clean.

- Connect the combination electrode to the meter.
- Measure with the combination electrode according to the operating manual of the meter and observe the following rules while doing so:

General rules for measuring

- TruLine ORP 51: Make sure the refilling opening for the reference electrolyte solution is open.
- Avoid displacement of the sample solution from one measurement to the next by taking the following measures:
 - Briefly rinse the calibration beakers with the solution the beakers are to be filled with next.
 - Between measurements, rinse the combination electrode with the solution specified below. Alternatively, you can also rinse the combination electrode with deionized water and then carefully dab it dry.
- Immerse the combination electrode in the solution in a vertical or slightly tilted position.
- Make sure the immersion depth is correct. The junction must be completely submerged in the solution. The junction is in the area of the bottom end of the shaft (see arrow).

The level of the reference electrolyte must be at least 2 cm above the level of the solution.

Conversion to the standard hydrogen electrode potential

$$U_H = U_{\text{Meas}} + U_{\text{Ref}}$$

where: U_H = ORP voltage, relative to the standard hydrogen electrode

U_{Meas} = measured ORP voltage

U_{Ref} = voltage of the reference system relative to the standard hydrogen electrode

U_{Ref} is temperature dependent and can be taken from the following table (also refer to DIN 38404-6):

T (°C)	U_{Ref} [mV] TruLine ORP 51	T (°C)	U_{Ref} [mV] TruLine ORP 51
0	+221	35	+187
5	+216	40	+181
10	+212	45	+176
15	+207	50	+171
20	+202	55	+165
25	+197	60	+160
30	+192		

Storage

During short measuring breaks

Immerse the electrode in the reference electrolyte.

Combination electrode	Reference electrolyte	YSI Model (see page 11)
TruLine ORP 51	3 mol/l KCl, Ag ⁺ -free	KCl-250 (250 ml)

Prior to the next measurement, briefly rinse the combination electrode with the test sample or deionized water.

Overnight or longer

Insert the clean electrode into the watering cap filled with reference electrolyte.

Note

During longer storage periods, salt deposits may develop on the watering cap. They do not affect the measuring characteristics and can easily be removed with deionized water when the combination electrode is put into operation again.

Aging

Every ORP combination electrode undergoes a natural aging process. Extreme operating conditions can considerably shorten the lifetime of the combination electrode. These are:

- Strong acids or lyes, hydrofluoric acid, organic solvents, oils, fats, bromides, sulfides, iodides, proteins
- High temperatures
- Great changes in pH and temperature.

The warranty does not cover failure caused by measuring conditions and mechanical damage.

Cleaning Remove water-soluble contamination by rinsing with deionized water. Remove other contamination as follows:

Contamination	Cleaning procedure
Fat and oil	Rinse with water containing household washing-up liquid
Lime and hydroxide deposits	Rinse with citric acid (10 % by weight)

After cleaning Rinse the combination electrode with deionized water.

Activating platinum electrodes


First-time activation during installation and as required

For the first-time activation use the activation powder from the SORT/RH reagent set (component of the SORT/RH reagent set). Immerse the moist (but not dripping) platinum electrode into the activation powder and rotate the electrode several times in the powder. Then remove the activation powder under flowing water with a soft brush (e.g. toothbrush).

Activation during very long set-up times

When changing from oxidizing to reducing test solutions and vice versa this can result in set-up times that can take significantly more than an hour. In this case pretreatment (activation) of the platinum surface can shorten the set-up time. The type of pretreatment (reducing or oxidizing) is based on the pH value and the ORP voltage (U_H) of the test solution where the latter must be estimated for the first measurement.

The type of pretreatment can then be determined using the following diagram where U_H is based on the normal hydrogen electrode:

Oxidizing pretreatment

Immerse the platinum electrode for two to three days in a sulfuric acid chlorine solution. Chlorine powder for producing the solution is included in the SORT/RH reagent set.

Note: The diaphragm must not be immersed in the chlorine solution!

Reducing pretreatment

When the electrode is ready for the test immerse it in the Zobell ORP buffer solution and wait for a stable measured value.

Wear parts and accessories

Description	Model	Order no.
Storage electrolyte solution 3 mol/l KCl, Ag ⁺ -free (250 ml)	KCl-250	109 705Y
Reagent set for regenerating ORP platinum electrodes, consisting of 10 g activation powder and 30 g chlorine powder	SORT/RH	109 730Y
Zobell ORP buffer solution (125 ml)	3682	061320

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Contact Information

Ordering & Technical Support

Telephone: (800) 897-4151
(937) 767-7241
Monday through Friday, 8:00 AM to 5:00 PM ET

Fax: (937) 767-1058

Email: environmental@ysi.com

Mail: YSI Incorporated
1725 Brannum Lane
Yellow Springs, OH 45387
USA

Internet: www.ysi.com

When placing an order please have the following information available:

YSI account number (if available)	Name and Phone Number
Model number or brief description	Billing and shipping address
Quantity	Purchase Order or Credit Card

Service Information

YSI has authorized service centers throughout the United States and Internationally. For the nearest service center information, please visit www.ysi.com and click 'Support' or contact YSI Technical Support directly at 800-897-4151.

When returning a product for service, include the Product Return form with cleaning certification. The form must be completely filled out for an YSI Service Center to accept the instrument for service. The Product Return form may be downloaded at www.ysi.com and clicking on the 'Support' tab.

What can Xylem do for you?

We're 12,500 people unified in a common purpose: creating innovative solutions to meet our world's water needs. Developing new technologies that will improve the way water is used, conserved, and re-used in the future is central to our work. We move, treat, analyze, and return water to the environment, and we help people use water efficiently, in their homes, buildings, factories and farms. In more than 150 countries, we have strong, long-standing relationships with customers who know us for our powerful combination of leading product brands and applications expertise, backed by a legacy of innovation.

For more information on how Xylem can help you, go to www.xyleminc.com

a xylem brand

YSI
1725 Brannum Lane
Yellow Springs, OH 45387
Tel: +1 937-767-7241; 800-765-4974
Fax: +1 937-767-1058
Email: environmental@ysi.com
Web: www.ysi.com

©Xylem Inc