

IMTEC-SALMONELLA-ANTIBODIES SCREEN (CUT-OFF)

Salmonella Screen

ELISA for the Detection of Antibodies against *S. typhimurium* and *S. enteritidis* (Ig(GAM))

Package Size

REF	ITC40040	96 Tests	Complete Testkit
IVD			

Please read the instructions carefully before testing.

Procedural precautions:

Do not use the reagents beyond the date of expiry.

[DIL] DB14, **[WASH]** 20x WB03, **[SUB]** TMB ELISA and **[STOP]** STOP ELISA may be interchanged between lots and test kits that share the same reagent designation.

All other reagents are specific for the individual test kit lot and must not be interchanged with other lots and test kits.

Store reagents at 2...8°C.

Intended Use

IMTEC-Salmonella-Antibodies Screen (cut-off) is an enzyme immunoassay (ELISA) for the qualitative measurement of IgG, IgM and IgA class autoantibodies against Salmonella in human serum. The assay is intended for in vitro diagnostic use only as an aid in the diagnosis of reactive arthritis.

The clinically most significant forms of reactive arthritis (ReA) develop within a few days to weeks after infection of the urogenital tract with chlamydia and the gastrointestinal tract with salmonella, yersinia, shigella or campylobacter respectively. Up to 40% of all cases, including a special form of Reiter's syndrome, are chronically progressive.

An early diagnosis based on a careful clinical case history and laboratory work up - including serological antibody detection - is therefore essential.

Up to date the Widal agglutination test is most frequently used for the diagnosis of salmonella infection, although it preferentially detects IgM antibodies typically appearing at the beginning of disease. The presence of high titres of IgA antibodies indicate the persistence of antigens in the intestine or joint and anti-IgG antibodies are reflecting a previous infection no longer active.

The IMTEC ELISA is detecting all three immunoglobulin classes IgG, IgM and IgA. Its sensitivity is superior to that of the Widal test, which detected only 38.5% of 130 patients with salmonella infections whereas 88.5% were detected by ELISA.

Principle

The test is based on the absorptive immobilisation of LPS from *Salmonella typhimurium* and *S. enteritidis* to the solid phase of microtiter strips and subsequent binding of anti-Salmonella antibodies from patient serum.

The bound antibodies are detected with a peroxidase-labelled secondary antibody that is directed against human IgG, IgM and IgA. After addition of substrate solution, a colour appears which intensity is proportional to the concentration and/or the avidity of the detected antibodies. Following the addition of stop solution, the colour switches from blue to yellow.

Reagents and Contents

[MTP]	12	Microtiter Strips (in 1 strip holder) 8-well snap-off strips, ready for use coated with LPS	
[CC]	2 ml	Cut-off-control (white cap), human, ready for use	
[NC]	2 ml	Negative Control Serum (green cap), human, ready for use	
[PC]	2 ml	Positive Control Serum (red cap), human, ready for use	
[WASH] 20x WB03	50 ml	Washing Buffer (black cap) Concentrate (20x) for 1 l TRIS buffer	pH 6.9 ± 0.2
[DIL] DB14	100 ml	Dilution Buffer (blue cap) ready for use Phosphate buffer	pH 7.3 ± 0.2
[CON]	15 ml	Conjugate Solution (white cap) anti-human-IgGAM HRP conjugate, ready for use	

[SUB] TMB ELISA	15 ml	TMB solution (black cap) ready for use, colourless to bluish 3,3', 5,5'-tetramethylbenzidin Hydrogen peroxide	pH 3.7 ± 0.2 1.2 mmol/l 3 mmol/l
[STOP] STOP ELISA	15 ml	Stop Solution (red cap) Sulphuric acid, ready for use	0.5 mol/l
	1	Adhesive Strip	

Safety Notes

Do not swallow the reagents. Avoid contact with eyes, skin and mucous membranes. All patient specimens and controls should be handled as potentially infectious. The controls have been checked on donor level for HCV and HIV-1/2 antibodies and HBsAg and found negative. Wear protective clothing and disposable gloves according to Good Laboratory Practices.

All materials contaminated with patient specimens or controls should be inactivated by validated procedures (autoclaving or chemical treatment) in accordance with applicable regulations.

Stability

The reagents are stable up to the stated expiry dates on the individual labels when stored at 2...8°C.

Reagent Preparation

Allow the testkit and all its components to reach room temperature before use! Used bottles should be closed carefully and stored at 2...8°C. Store **[SUB]** protected from light.

Do not use polystyrene vessels for handling of **[CON]**.

To avoid potential microbial and/or chemical contamination, unused reagents should never be transferred into the original vials.

Washing Buffer Solution **[WASH]**

Any crystallised salt inside the bottle must be resolved before use. Dilute 1 part **[WASH]** 20x with 19 parts distilled water. **[WASH]** is stable for 6 weeks stored at 2...8°C.

Specimen

Patient sera

Use samples freshly collected or freeze samples at -20°C. **Freeze and thaw once only.** Do not use serum samples inactivated by heat treatment at 56°C.

Allow the samples to reach room temperature (30 min.).

Dilute sera 1:101 with **[DIL]** (add 10 µl serum to 1 ml **[DIL]**).

Procedure

- **Pipette 100 µl** diluted sample, **[CC]**, **[PC]** and **[NC]** into **[MTP]**, for blank use **[DIL]** instead of sample dilution, seal **[MTP]** with adhesive strip.
- Incubate for **1 hour** at RT.
- Discard the solution from **[MTP]**. Wash **[MTP]** 3 times using 300 µl **[WASH]** per well.
- Discard **[WASH]** and knock out residues on an absorbent paper or cloth.
- **Pipette 100 µl** **[CON]** and seal **[MTP]** with adhesive strip.
- Incubate for **30 min.** at RT.
- Discard the solution from **[MTP]**. Wash **[MTP]** 3 times using 300 µl **[WASH]** per well.
- Discard **[WASH]** and knock out residues on an absorbent paper or cloth.
- **Pipette 100 µl** **[SUB]** and incubate for **10 min.** At room temperatures above 25°C the substrate incubation could be shortened, but should never fall short of 5 min..
- **Add 100 µl** **[STOP]** per well.
- Read absorbance values at 450 nm within the next 10 min. after stopping. Bi-chromatic measurement with a reference wavelength at 620 – 690 nm is recommended.

Automation

The IMTEC-Salmonella-Antibodies Screen (cut-off) ELISA may be processed with suitable automated ELISA analyzers. Applications have to be validated prior to diagnostic use.

Validation of the Test

The test results are valid provided that the following criteria are met for the obtained results:

- $PC > CC$.
- PC does not fall below an absorbance value of 0.4.
- $NC < CC$.
- $PC / CC = 1.2 - 5.0$

In order to improve accuracy of the test results we recommend to run PC , CC , NC and patient samples in duplicate.

Interpretation of results

Interpretation of results can be made by comparing the absorbances of CC and of the samples:

- Absorbances $> 1.1 \times CC$ have to be considered as positive.
- Absorbances $< 0.9 \times CC$ have to be considered as negative.
- Absorbances $\geq 0.9 \times CC$ and $\leq 1.1 \times CC$ have to be considered as equivocal.

Limitations

A positive result must be used in association with clinical evaluation and diagnostic procedures. The values obtained from this assay are intended to be an aid for diagnosis only.

Elevated anti-Salmonella antibodies may occur in individuals with no evidence of clinical disease.

If the patient sample contains elevated levels of immune complexes or other immunoglobulin aggregates, false positive results by non-specific binding cannot be ruled out.

The performance characteristics for this assay have not been established for plasma samples.

Performance Characteristics

Typical performance data can be found in the Verification Report, accessible via:

www.human.de/data/gb/vr/el-40040.pdf or

www.human-de.com/data/gb/vr/el-40040.pdf

If the performance data are not accessible via internet, they can be obtained free of charge from your local distributor.

Safety Notes

STOP Warning

• Hazard statements

H315 Causes skin irritation.

H319 Causes serious eye irritation.

SUB Danger

• Hazard statements

H360D May damage the unborn child.

• Precautionary statements

CC NC PC WASH 20x DIL CON SUB STOP

P234 Keep only in original container.

P260 Do not breathe dust/fume/gas/mist/vapours/spray.

P262 Do not get in eyes, on skin, or on clothing.

P281 Use personal protective equipment as required.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337+P313 If eye irritation persists: Get medical advice/attention.

P401 Store in accordance with local/regional/national/international regulations.

P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

References

1. Isomäki O. *et al.*, Lancet June **24**, 1411-1414 (1989)
2. Mäki-Ikola O. *et al.*, J. Infect. Dis. **164**, 1141-1148 (1991)

EL-40040 INF ITC40040 GB 03-2018-13



Human