

Digital Refractometer for Ethyl Alcohol PET-109



The concentration of ethyl alcohol used for foods, cleaning, cooling and other purposes can be measured with only a small quantity of sample and with a simple procedure. The temperature is automatically corrected and not necessary to worry about the temperature.

Specifications

Model	PET-109
Cat.No.	3486
Scale	Concentration of ethyl alcohol
Range	0.0 to 45.0% (W/W)
Resolution	0.1%
Accuracy	0 to 30%:±0.5% Around 40%:±0.8% Around 45%:±1.2%
Temperature compensation	10 to 35°C
Power Supply	006P dry battery (9V)
Dimensions & Weight	17×9×4cm, 270g

Options

* Calibration Certificate : Contact an ATAGO representative for details.