

Metallographic Sample Grinding & Polishing Machine

Model: JG-162



(For reference only)

Features and Application

- Suitable for metallographic specimen pre-grinding, grinding and polishing operation.
- Equipped with cooling device to prevent the sample from overheating
- Left disc used for grinding and right disc used for polishing, so two people can operate it at the same time.
- Easy to operate, safe and reliable. It is an ideal sample preparation equipment for factories, research institutes and universities laboratory.



Technical Specifications

Model	JG-162
Grinding disc diameter	230mm
Polishing disc diameter	200mm
Grinding speed	450r/min
Polishing speed	600r/min
Motor driver	370W, 220V, 50Hz
Dimension (L x W x H)	700 x 600 x 278mm
Weight	45Kg

Details

