


### NEXUS WIRELESS CALL SYSTEM



The spread of coronavirus has led to a significant surge in demand for medical ward services. In order to cope with bed short of the expected surge in corona patients, hospitals have already begun to step up and provide overflow bed capacity, and the call bell is an indispensable part of it. The *Nexus Call System* from Rondish may come in handy. It is a portable nurse call solution that *can be quickly deployed with mimnimum installation*. It's high flexibility making it a perfect choice for a surge in bedding needs during the coronavirus pandemic.









#### SCHEMATIC DIAGRAM

#### SCHEMATIC DIAGRAM - RONDISH NEXUS WIRELESS NURSE CALL SYSTEM



SYSTEM LAYOUT



**Eightfold Sdn Bhd** (775128-W) 69-1, Jalan SS19/6, 47500 Subang Jaya, Selangor Darul Ehsan, Malaysia Tel: 03-5612 8169 www.eightfold.com.my



### EXPANDABILITY

### Fall Management Components



Wired Bed Monitoring





### **Cordless Bed Monitoring**





#### **TXP-11 Pendant Transmitter**

- Active emergency button sends signal to central monitor or field receiver
- Acceleromater to trigger call in case of patient fall



**PRODUCT DATASHEET** 

#### Product

#### : NDP-11 Nexus Display Panel

- Description
- 2 Android tablet with Nexus software for basic call system. Must include receiver dongle to operate

#### Features:

- 10" Android tablet with folio included (as shown) •
- Customizable descriptions for device location •
- Three alarm priority levels follows HTM08 •
- Built-in record of calls
- Low battery and signal lost alert
- 6 hours battery power backup in case of main • failure

433.92MHz

5°C to 40°C

850g

100m (line of sight)

266 x 155 x 20mm

DC5V input, 2000mA Battery back-up for 4+ hours

Wall-mount bracket available

#### Works with:

- DON-20 Wireless Receiver Dongle
- ZP-01(M) Master Protektor Unit
- **NEC-05** Wireless Patient Handset
- WCP-11 Wireless Call/Reset Button

#### **Technical Specifications**

Operating frequency:

Operating range:

Operating temp.:

Power supply:

Dimensions:

Net weight:



- WCP-11C + NEC-12 Wireless Call/Reset with Patient Handset
- ILB-11 Room Indicator Light
- TXP-11 Emergency Transmitter with Man-Down alert
- WLSMCL-01 Wireless Multi-colored Light

#### Standards

EMC:	EN 55022
Radio:	EN 300220-1&2,
	EN 301489-1&3
Alarm:	EN 50130-4
Safety:	EN 60950-1
Manufacture:	ISO 9001:2008, ISO 13485:2003

#### Product

Description

#### NEC-02 / NEC-12 Nurse Call Handset 2

Nurse call handset with call button and wall-mount bracket. Connects with Multifunction or Bedhead Call Point

#### Features:

- Large rubber call button with reassurance LED
- USB break point to protect products
- NEC-02 cable comes with RJ12 connection for Protektor Call Point. NEC-12 comes with RJ10 connection instead
- Plastic wall-mount bracket with optional 90° 6.3mm jack lead

#### Works with:

- WLSCP-02 Multi-function Call Point (NEC-02)
- WCP-21C Bedhead Call Point (NEC-12)
- MEX-20 Multiport Expander
- CMEX-11 Cordless Expander
- USBC-05 Lead with 90° audio jack (6.3mm)

### **Technical Specifications**

Dimensions(LxWxD) ~12 x 5 x 2.5cm (handset only) Housing material: ABS 160g (with cable and bracket) Net weight: 5°C to 40°C Operating temp.: Rel. Humidity: 85%



2

2

#### Eightfold Sdn Bhd (775128-W) 69-1, Jalan SS19/6, 47500 Subang Jaya, Selangor Darul Ehsan, Malaysia Tel: 03-5612 8169 www.eightfold.com.my



#### Product

#### WCP-21C Bedhead Call Point

Description

Wireless call point with both call and reset buttons, with RJ10 input for a patient handset, works as part of Protektor network

#### Features:

- Handshake signaling with Protektor / Nexus network
- Read-in to indication light
- Generate call, reset, emergency call and nurse presence
- Optional patient handset (as shown)
- Programmable to display variety of locations
- Low battery and automatic fault detection
- Changeable batteries
- Battery life: ~2 3 year, assumes 3-5 calls/day

#### Works with:

- Nexus Display Panel
- WLSMCL-01 Wireless Multicolored Light
- ZP-01 Zone Processor
- NEC-12 Patient Handset (shown in picture)

#### Product : WCP-21 Bathroom Call Point

Description

Wireless call point with call and reset buttons, waterproof to IP 65, and works as : part of a Protektor network

#### Features:

- Handshake signaling with Protektor network
- Read-in to indication light
- · Safe for use in toilet or shower area
- Normal call (red button) and Reset (green button)
- Optional pull cord (not shown)
- · Programmable to display variety of locations
- Waterproof to IP65
- Low battery and fault detection
- Changeable batteries
- · Battery life: 2-4 years, assumes 1 call/day

#### Works with:

- NDP-11 Nexus Display Panel
- WLSMCL-01 Wireless Multicolored Light
- ZP-01 Zone Protektor

#### **Technical Specifications**

Operating frequency: 866 or 922MHz			
Operating range:	200m (line of sight)		
Power supply:	2 x CR2477 lithium ion battery		
Dimensions:	8.5 x 8.5 x 6cm		
Housing material:	ABS		
Net weight:	130g		
Operating temp .:	5°C to 40°C		
Rel. humidity:	95%		



EMC:	EN 55022
Radio:	EN 300220-1&2,
	EN 301489-1&3
Alarm:	EN 50130-4
Safety:	EN 60950-1
Manufacture:	ISO 9001:2008, ISO 13485:2003





#### Product

#### : WLSCP-02 Multi-function Call Point

Description

: Wireless multi-function call point with five buttons, input for patient handset, and bed sensorpad/floor mat

#### Features:

- Handshake signaling with Protektor network
- Normal call (red button)
- Emergency call (yellow button)
- Presence1/Presence2 (grey button)
- Guard presence (blue button)
- Reset (green button)
- Input for nurse call handset, RJ12 connection
- Input for bed pad/floor mat, RJ10 connection
- Comes with removable pull cord
- Low battery and fault detection

#### Works with:

- WLSMCL-01 Wireless Multicolored Light
- ZP-01 Zone Protektor
- NEC-02 Nurse Call Handset
- FBP-10 Foldable Bed Pad
- FMAT-02 Floor Mat
- CK-01 Caregiver Key
- WPG-01 System Programmer
- NDP-11 Nexus Display Panel

#### **Technical Specifications**

Operating frequency: 866 or 922MHz			
Operating range:	200m (line of sight)		
Power supply:	12VDC, 200mA		
	4 x AA alkaline batteries		
Dimensions:	12.5 x 10 x 3.5cm		
Housing material:	ABS + PC		
Net weight:	200g		
Operating temp.:	5°C to 40°C		
Rel. humidity:	< 90%		



#### Standards

EMC:	EN 55022
Radio:	EN 300220-1&2,
	EN 301489-1&3
Alarm:	EN 50130-4
Safety:	EN 60950-1
Manufacture:	ISO 9001:2008, ISO 13485:2003
FCC Approved	FCCID: WNG-WLSCP-02



#### Product : WLSMCL-01 Wireless Multi-coloured Light

Description

: Wireless indication light with four colors to indicate call type, can forward signals received to boost signals and increase range, optional buzzer

#### Features:

- Four colors, corresponding to normal call, toilet alarm, and two levels of nurse presence
- Optional buzzer setting for active call
- Red or white light flashes for emergency call
- Wireless receiver and transmitter allows for signal- boosting to increase wireless range
- Flashing blue LED to indicate power status

#### Works with:

- WLSCP-02 Multi-function Call Point
- WCP-21C Bedhead Call Point
- WCP-21 Waterproof Call Point
- ZP-01 Zone Protektor

#### **Technical Specifications**

Operating frequency: 866 or 922MHz Operating range: 100m (line of sight) Power supply: 12Vdc, 200mA Battery life: 2-4 years Dimensions: 16 x 12 x 4cm Housing material: ABS + PC Net weight: 180g Operating temp.: 5°C to 40°C Rel. humidity: <90%



#### Standards

EMC:	EN 55022
Radio:	EN 300220-1&2,
	EN 301489-1&3
Alarm:	EN 50130-4
Safety:	EN 60950-1
Manufacture:	ISO 9001:2008, ISO 13485:2003
IP Rating	65
FCC Approved	FCCID: WNG-WLSMCL-01



#### Product **GP2009TR POCSAG Transmitter / Repeater** :

Description

POCSAG transmitter/repeater. Sends alarm details to wireless LED display : panel, dotmatrix display and/or data pager.

#### Features:

- Transmits alarm data to remote displays
- Simple software interface
- Adjustable data rate on transmission
- Can repeat signals to extend site coverage
- Wall or table mountable, adjustable antenna
- USB and RS232 ports

#### Works with:

- GP2015S wireless LED display panel
- GP2009N data pager •
- Protektor II Wireless System
- Nexus Wireless Call System
- SmartBeds Ward Management Platform

#### **Technical Specifications**

Operating frequency: 420 or 470MHz Operating range: Up to 2km Transmission Power: 2W max Power supply: 7-20Vdc Dimensions: 110 x 62 x 22mm



Standards		
EMC:		
Radio:		
Alarm:		

RoHS FCC

CE

Housing material:	Aluminum
Net weight:	230g
Operating temp.:	-10°C to 50°C
Rel. humidity:	<90%

#### Product : GP2015S Wireless LED Dotmatrix Display

Description Lightweight wireless LED dot matrix display for showing alarm type and location. : Works with Protektor II, Nexus and SmartBeds

#### Features:

- 2 x 16-digit dot matrix alphanumeric display
- Wall or ceiling mounted •
- Audible buzzer or voice alert
- Supports English, German, French, and Chinese
- 8 capcodes allows multiple zones in a facility •

#### Works with:

- **GP2009TR** transmitter
- Protektor II Wireless System
- Nexus Wireless Call System

SmartBeds Ward Management Platform

Technical Specifications					
Operating frequency:	400 - 499MHz	Housing material:	ABS	EMC:	CE
Operating range:	Up to 2km	Net weight:	1600g	Radio:	RoHS
Transmission Power:	N/A	Operating temp.:	-5°C to 40°C	Alarm:	FCC
Power supply:	5Vdc, 3A	Rel. humidity:	<90%		
Dimensions:	650 x 168 x 32mm				

