

## Denture Base Polymers

Through particular multidirectional rotation compounding process, this product is made up by powder material of imported high-quality polymethyl methacrylate resin, with cadmium-free inorganic red pigment and bionic stria vascularis added.

This product is clinically suitable for prosthodontics as denture base of patients with hypodontia.

### Main Feature

- Easy to polish
- Excellent handling properties
- Blanch Free
- Craze free and Solvent resistant
- Good working time
- Cadmium free



### Classification

Type	Mixing Ratio	Mixing Time	Dough Time	Working Time	Curing Time
Fast Heat Curing	1ml Liquid:2.4g Powder	1-2mins	14-20mins	8-10mins	20mins
Traditional Heat Curing	1ml Liquid:2.34g Powder	1-2mins	15-20mins	8-10mins	90mins
Self Curing	1ml Liquid:2.2g Powder	0.5mins	12-14mins	5-6mins	30mins

### Package

Powder: ◦ 100g/bag, 10 bags /box;

◦ 1000g/can;

Liquid: ◦ 500ml/bottle/box



### Shade

- 5 Shades: 2S,2ST,3ST, I, V
- Popular Shades: 2S, 2ST,3ST
- Customized shade also available

