

Beta-Human Chorionic Gonadotropin (Beta-HCG) Test Kit (Colloidal Gold)



Applications:

Professional institutions, hospitals, clinics, drugstore and personal home usage, etc.

Technical Parameters:

Methodology	Colloidal gold method
Sample Type	Urine
Storage Condition	4-30°C
Period of Validity	24 months
Sample Volume	3 drops, about 80µL
Detection Duration	5-10 minutes
Detection Limit	10mIU/mL
Package Information	20pcs/box, 30boxes/carton
Packing Size	310*330*360mm
Gross Weight	8kg