

## Typhoid Rapid Test Kit (Colloidal Gold)



### Applications:

Hospitals, customs, health and epidemic prevention departments, drugstore and personal home usage, etc.

### Technical Parameters:

Methodology	Colloidal gold method
Sample Type	Whole blood/serum/plasma
Storage Condition	4-30°C
Period of Validity	24 months
Sample Volume	20µL Whole blood/serum/plasma
Detection Duration	15-20 minutes
Performance Index	The sensitivity was 83.5%, the specificity was 95.2%, and the total coincidence rate was 89.9%.
Package Information	20pcs/box, 30boxes/carton
Packing Size	310*330*360mm
Gross Weight	8kg