

Chikungunya IgG/IgM Combo Rapid Test Kit (Colloidal Gold)



Applications:

Hospitals, customs, health and epidemic prevention departments, drugstore and personal home usage, etc.

Technical Parameters:

| Methodology | Colloidal gold method |
|---------------------|---|
| Sample Type | Whole blood/serum/plasma |
| Storage Condition | 4-30°C |
| Period of Validity | 24 months |
| Sample Volume | 20µL Whole blood/serum/plasma |
| Detection Duration | 15-20 minutes |
| Performance Index | The sensitivity was 88.2%, the specificity was 96.3%, and the total coincidence rate was 91.8%. |
| Package Information | 20pcs/box, 30boxes/carton |
| Packing Size | 310*330*360mm |
| Gross Weight | 8kg |