

# BOT

## BRITISH OXYGEN TECHNOLOGY



**BOT**  
BRITISH OXYGEN  
TECHNOLOGY



**5L Medical  
Oxygen  
Concentrator  
FOLEE**

# Y007C-5W



**011-64806728 ELSON**

# BOT BRITISH OXYGEN TECHNOLOGY



OUR

HOUSEHOLD OXYGEN GENERATOR

## 5L MEDICAL GRADE OXYGEN-CONCENTRATOR

Real-time display of oxygen flow and concentration

Foylee 富林  
Y007C-5W

- O<sub>2</sub> Function of releasing anion
- Easy to use
- Mute design
- Oxygen production and atomization
- Multiple purification

REAL-TIME MONITORING, SAFETY PROTECTION

## LOW OXYGEN / FAULT ALARM

1. OXYGEN PURITY LESS THAN 82%, YELLOW LIGHT ALARM
2. EQUILIBRIUM SYSTEM WARNING WHEN OXYGEN SYSTEM IS ABNORMAL.
3. RED LIGHT ALARMS WHEN PRESSURE IS ABNORMAL
4. THE MACHINE BEEPS WHEN THE POWER IS ABNORMAL

LOW OXYGEN

PRESSURE CYCLE FAULT

OUT OF POWER

COMPRESSOR FAULT

STRENGTH

5L Medical Oxygen Concentrator Y007C-5W



011-64806728 ELSON









Distributed by : British Oxygen Technology | 3, Jalan Permatang 8, 81100 Johor Bahru, Johor

# BOT

# BRITISH OXYGEN TECHNOLOGY



FUNCTIONS & FEATURES

 NEBULIZER	 REMOTE CONTROL
 OXYGEN FLOW FINE ADJUST	 TIMING FUNCTION
 93% OVER 93% HIGH OXYGEN PURITY	 MUTE DESIGN
 LED DISPLAY	 FAULT ALARM

## 8 USEFUL ADVANCED FEATURES

OXYGEN & NEBULIZER **ONE MACHINE DUAL USE**

OXYGEN PURITY OVER 93% AT ALL THE FLOW

## OUR

## STRENGTH

5L Medical Oxygen Concentrator Y007C-5W



011-64806728 ELSON

Distributed by : British Oxygen Technology | 3, Jalan Permatang 8, 81100 Johor Bahru, Johor

# BOT

# BRITISH OXYGEN TECHNOLOGY



## 5L MEDICAL OXYGEN CONCENTRATOR

### PARAMETER

MODEL NUMBER	Y007C-5W
OXYGEN PURITY	93% +/- 3%
FLOW	0L - 5L/MIN
POWER	300W
OUTPUT PRESSURE	0.02 - 0.07MPA
VOLTAGE	AC220V + 22V, 50Hz + 1Hz
DIMENSION	400*260*595MM
PACKAGING	430-380-67.5MM
NET WEIGHT	15KG
TOTAL WEIGHT	18KG
NOICE	<60dB(A)
ANION	5000000/CM3
DISPLAY CONTROL	LED DISPLAY
TECHNOLOGY	PSA

# DATASHEET

5L Medical Oxygen Concentrator **Y007C-5W**



**011-64806728 ELSON**

Distributed by : British Oxygen Technology | 3, Jalan Permatang 8, 81100 Johor Bahru, Johor

# BOT BRITISH OXYGEN TECHNOLOGY

**MDA APPROVED**



LAMPIRAN 1  
Attachment 1



No. Pendaftaran: **GB8276720-39096**  
Registration No.:

Butir-butir peranti perubatan yang didaftarkan  
Particulars of the registered medical device

Nama Peranti Perubatan **OXYGEN CONCENTRATOR**  
Medical Device Name

Kelas **CLASS B** Jenama **Y007 SERIES**  
Class Brand

Kelompok **SINGLE** Produk Identifikasi **Y007C-5W**  
Group Product Identifier

Nama dan alamat pembuat: **JIANGSU FOLEE MEDICAL EQUIPMENT CO., LTD.**  
Name and address of manufacturer **NO. 16 XINGMAO ROAD, 212009 ZHENJIANG CITY, JIANGSU PROVINCE, CHINA, 212009 CHINA**

Medical Device  
AUTHORITY  
MALAYSIA

# Tips of choose POC

---

## How to choose a suitable device for your family member?

### 1. What oxygen should be Enough for your family member?

Usually doctor will write or print oxygen dose and oxygen therapy time in the prescription. If doctor write 2LPM or 2 liters/min in prescription, that means 2 liters oxygen with above 90% oxygen purity.

If doctor write 3LPM or 3 liters, that means 3 liters oxygen with above 90% oxygen purity.



### **Be Attention:**

Please make confirmation that machine could produce 90% saturation oxygen at least. There are many devices just produce 70% saturation oxygen or less when you turn up flow. ( oxygen saturation less than 90% is not allowed by medical oxygen therapy. Insufficient oxygen will worsen illness and be harmful to your family member life)

Another thing to consider is choose bigger oxygen flow than your family member actually need. If his oxygen prescription is 3LPM, we suggest you choose a device which has capacity of above 3 liters flow with 90% oxygen saturation. Your family member might need more oxygen in winter and in future. You need to take it into consideration.

## 2. Put safety of your family member first.

Low oxygen alarm could monitor oxygen saturation. If machine deliver oxygen less than 82% saturation, alarm will work to remind you. If alarm activate, you need to call Lovego service station to check machine.



### Choose a device with low oxygen alarm:

If you choose a device which doesn't have low oxygen alarm, you won't know whether it is producing oxygen or normal air. If device has problem in oxygen concentration, it will just deliver normal air. But device still looks working normally. You have no idea when your family member inhale normal air. Your family member will feel suffocating. It will worsen illness and be harmful to your family member life. It will cause syncope and death in emergency situation.

A man in a blue suit with arms crossed, standing in front of a blurred background.

# BRITISH OXYGEN TECHNOLOGY

TOTAL GAS SOLUTION

**011-64806728 Elson**