

TO OPEN LID,
PRESS
WHITE
BUTTON



RESCUE
READY®


Defibrillation Electrodes
Electrodes de défibrillation
Elettrodi di defibrillazione
Defibrillator-Elektroden
Electrodos
Electrodos
自動除去用電極
自動電極

REF 9131

POWERHEART® AED G3

automated external defibrillator

XXXX-XX-XX
LOT 170327-01
2017-03-27

 **CARDIAC
science**

Manufactured in USA for:
Cardiac Science Co.
N7 W28025 Johnson Dr.
Waukegan, WI 53186 U.S.
www.cardiacscience.com

MDSS GmbH
Schiffgraben 41
D-30175 Hannover
Germany

 **CARDIAC
science™**

U.S. Patents
5,402,884; 5,579,819;
5,984,102; 4,926,379

R ONLY CE G006

TEAR ALONG DOTTED LINE

TEAR HERE TO OPEN

Powerheart® G3 Elite AED

KEY FEATURES

EFFECTIVENESS

- RHYTHMx® - a proven technology which analyses the patient's heart rhythm in seconds
- STAR® biphasic technology - delivers a customised shock at an appropriate energy level for each patient
- Paediatric capability

SIMPLICITY

- RescueCoach™ - provides user-paced step-by-step voice and text prompts
- Bright backlit text display - displays rescue prompts and critical rescue information
- Non-polarised defibrillation electrodes for quick placement on either side of the patient's chest
- CPR - simple coaching, rate metronome

CONFIDENCE

- Rescue Ready® technology – self tests device, pads, and battery for presence and functionality every day
- Rescue Ready® indicator
- 7-year device warranty
- 4-year full operational performance guarantee on battery*
- Strong indemnification policy

* Some conditions apply. Contact your Cardiac Science representative for details.

WHY WOULD YOU NEED AN AED?

To treat a Sudden Cardiac Arrest (SCA) - for when a person is unresponsive and not breathing or not breathing normally.

TARGET SECTORS

Education
Government
Emergency Services

PRODUCT CODES

9790A-1002 - Fully Automatic
9790E-1002 - Semi-Automatic
9730-002 - Paediatric Defibrillation Electrodes

COMPETITIVE DEVICES INCLUDE

Heartsine 350P & 360P | Mediana A15 | CU Medical iPad



If you have any questions or queries about any of the above please contact Cardiac Science
Tel: +44 (0)161 926 0000 • Fax: +44 (0)161 926 0002
UK Sales: sales@cardiacscience.co.uk
UK Technical Support: uktechsupport@cardiacscience.com

MKT-00273-02rA