


MIX Medical Cart [Procedure]


The **MIX PROCEDURE CART** can be used in multitude of medical departments such as:

1. supply/storage
2. Med/Surgical
3. Cast
4. Dressing
5. Airway
6. Bio-Medical
7. I.V. Therapy
8. Critical Care
9. Trauma
10. Lab



IV Pole

IV Pole can be adjusted




Cart size options



Cart specifications

Item No	W		D		H		w/o Casters	
	mm	inch	mm	inch	mm	inch	mm	inch
MX37PRO	690	27-1/8"	518	20-3/8"	1085	42-11/16"	951	37-7/16"
MX34PRO	690	27-1/8"	518	20-3/8"	1006	39-5/8"	872	34-5/16"
MX31PRO	690	27-1/8"	518	20-3/8"	927	36-2/1"	793	31-3/16"
MX28PRO	690	27-1/8"	518	20-3/8"	848	33-8/3"	714	28-1/8"




MX34PRO



Different accessory combinations

Accessory can be set up and adjusted easily and freely by requirements


Accessory package (standard)



IV Pole



Waste container set(w/ lid)



Sharps container set



Multi storage (L)