


Linovera[®]

Solution of hyperoxygenated fatty acids

B. Braun Wound Care. Right. From the start.



Wound Management

Linovera®

Solution of hyperoxygenated fatty acids

Protects against external agents

Improves skin resistance

Repairs the upper epidermis

Restores capillary circulation

Stimulates the synthesis of collagen

Has a hydrating and healing effect



Linovera® is a solution of hyperoxygenated fatty acids used in the prevention and treatment of Stage I pressure ulcers. Due to the high percentage of linoleic acid, its repairing action is increased.

The vegetable extracts – aloe vera and centella asiatica – have a major protecting and healing effect on the skin.

Composition

- › Hyperoxygenated essential fatty acids (linoleic acid)
- › Tocopherols
- › Aloe Vera
- › Centella Asiatica
- › Perfume



Indications

- › Prevention and treatment of Stage I pressure ulcers

TIP FOR APPLICATION

- › Spray Linovera® 2 – 3 times a day on the affected area and rub in softly to evenly distribute the product until it is fully absorbed
- › Use in combination with Askina® Heel for optimal protection of the heel area against pressure ulcers

Linovera®	Size	Pcs/Pack	Article number
GB ES PT IT	30 ml	1	467933
DE CZ PL SK			468156
FR NL RU DK			468157
FR DE NL IT			468155
GB ES PT IT	30 ml	10	468160
DE CZ PL SK			468161
FR NL RU DK			468166
GB ES PT IT	50 ml	10	468158
DE CZ PL SK			468159
FR NL RU DK			468162