

Metallographic Sample Polishing Machine

Model: JG-153



(For reference only)



Features and Application

JG-153 Metallographic Specimen Polishing Machine is suitable for the polishing of the specimen which has been grinded. The specimen surface is very smooth after processing and can be used by people to observe and measure the metallographic structure of specimen under microscope.

Technical Specifications

Model	JG-153
Diameter of polishing disc.	200mm
Polishing disc rotating speed	1400rpm
Motor power	180W
Power supply	220V, 50Hz
Dimensions	45 x 39 x 28cm
Net weight	16Kg

