

# Metallographic Sample Grinding & Polishing Machine

Model: JG-161



(For reference only)

## Features and Application

- Double-disc desktop, can be operated by two people simultaneously;
- Regulating speed by a frequency changer, with a speed of 50 ~ 1000rpm;
- Equipped with a cooling device, preventing the damage of the metallographic structure caused by overheat;
- Applicable to pre-grinding, grinding and polishing of metallographic samples;
- Easy to operate, safe and reliable, is an ideal equipment for the labs of plants, research institutes and universities and college.



## Technical Specifications

Model	JG-161
Diameter of grinding disc	200mm (203mm is optional)
Grinding disc rotating speed	50 ~ 1000rpm
Diameter of polishing disc	200mm
Polishing disc rotating speed	50 ~ 1000rpm
Working voltage	220V / 50Hz
Diameter of Abrasive Paper	200mm
Motor	YSS7124, 550W
Dimension	700 x 600 x 278mm
Weight	50Kg

## Details

1. With cabinet (optional)



## 2. Panel

