Paperboard Compression Tester

Model: JG-300



(For reference only)

Equipment Usage

The instrument is mainly used for: ECT, RCT, PAT, FCT, honeycomb cardboard, etc...

Standard

• GB/T 6546, GB/T 6548, GB/T 2679.8, GB/T 2679.6, GB/T 22874, GB/T 6544, ISO 3035, ISO 3037, ISO 7263, TAPPI-T472, T808, T810, T811, T818, JIS-P8126 etc





Features

- Product using microcomputer automatic control, easy operation, high accuracy, stable and reliable performance.
- Products using NXP (PHILPS) high performance ARM core processor, the United States Logic Cirrus dual channel 24 bit independent high-speed sampling chip, Japan OMRON electronic components, to ensure accurate sampling and high-speed operation, to ensure the accuracy of the whole.
- Products using 5.1 inch LCD screen, the English menu operation interface, test: real-time display of the test value, force value change curve and other information.
- Automatic measurement sample, automatic lock peak, powerful data display and analysis management capabilities, 99 sets of data storage, power down or power data memory function, query function, print function, stroke protection, overload protection, automatic reset, fault alarm, sensor nonlinear correction and other functions, effectively ensure the safety and simplicity of user operation.
- N, kgf, lbf switching between multiple units of measure, conversion function.
- Products with high speed micro printer, print content are: test time, test number, test data and related standards set of test indicators.
- Real time clock, replacing the previous manual input date or manual fill in the date of the tedious operation.

Technical Specifications

Model	JG-300	
Range	3000N (300Kg)	
Error of indication value	<1%	
Resolution	0.1N or 1N	
Test speed	12.5 ± 2.5mm/min	
Pressing plate diameter	120mm	
Maximum stroke	80mm	
Power supply	AC220V 3A 50Hz	
Dimensions	360mm(L) x 330mm(W) x 510mm(H)	
Weight	55Kg	

According to the user's own needs, choose and purchase related accessories:

- 1. ECT: JG-349 ECT(PAT)Sample Cutter
- 2. RCT: JG-350 RCT Sample Cutter, JG-355 RCT Accessories
- 3. PAT: JG-349 ECT(PAT)Sample Cutter, JG-356 PAT Accessories
- 4. FCT: JG-360/JG-361 FCT Sample Cutter



Optional

JG-349 ECT (PAT) Sample Cutter







Equipment Usage

• Products are mainly used for cutting ECT, PAT specimen.

Standard

GB/T 6546, GB/T6548, GB/T 6544, QB/T 1671 etc.....

Features

- The thrust type sampling, the structure is simple and firm.
- The double stainless steel ball bearing guide sleeve design, equipped with imported Japanese blade demolition convenient, sampling effort, accurate, without burr.

Model	JG-349 ECT(PAT) Sample Cutter	
Maximum sample thickness	18mm	
Sample size	100mm, 80mm, 25mm (ECT sample 100 x 25mm; PAT sample 80 x 25)	
Sampling accuracy error	±0.5mm	
Dimensions	440mm(L) x 380mm (W) x 180mm(H)	
Weight	19Kg	



JG-350 RCT Sample Cutter





Equipment Usage

Products are mainly used for cutting RCT, CCT, CMT specimen.

Standard

• GB/T 2679.8, QB/T 1671 etc.....

Features

- The use of stampling, simple structure, strong.
- Double ball guide sleeve design, blade after quenching, sampling effort, accurate, no burr.

Model	JG-350 RCT Sample	
Sample thickness	0.1mm ~ 1.0mm	
Sample size	152±0.2mm x 12.7±0.1mm	
Dimensions	530mm(L) x 150mm (W) x 405mm(H)	
Weight	35Kg	






JG-355 RCT Accessories







Equipment Usage

• Product is special attachment to RCT.

Standard

• GB/T 2679.8 etc......

Features

- Products made of stainless steel material, fine workmanship, the use of light, convenient.
- Using a wide range, a total of 13 central disk, specifications 0.15 ~ 1.00mm, with a different center plate, can be adapted to different thickness of the specimen test.

Model	JG-355 RCT Accessories	
Central disk	13 central disk, specifications 0.15 ~ 1.00mm	
Outer packing box size	330mm(L) x 230mm (W) x 50mm(H)	





JG-356 PAT Accessories





Equipment Usage

• Product is special attachment to PAT.

Standard

• GB/T 6548, GB/T 6544 etc.....

Features

- Fine workmanship, the use of light, convenient.
- The use of a wide range, can be used with different corrugated paper test.

Model	JG-356 PAT Accessories	
Applicable type	A-type, B-type, C-type	
Pin	2mm, 3mm	
Outer packing box size	200mm(L) x 190mm(W) x 50mm(H)	
Weight	2Kg	



JG-360/JG-361 FCT Sample Cutter





Equipment Usage

• The product is mainly used for cutting the FCT sample, and the sample size is divided into two sections. Two kinds of specification choose one.

Standard

• GB/T 2679.6, QB/T 1671 etc.....

Features

- Compact structure, easy to carry.
- By rotary sampling, convenient sampling, no burr, the blade can be easily removed and replaced (with spare blade gift).

Model	JG-360 FCT Sample Cutter	JG-361 FCT Sample Cutter
Sampling area	32.2cm² (64mm±0.2mm)	64.5cm ² (90.6mm±0.5mm)
Maximum sample thickness	<15mm	<15mm
Dimensions	130mm(L) x 130mm(W) x 170mm(H)	150mm(L) x 150mm(W) x 170mm(H)
Weight	2.5Kg	3Kg

