


JG-349 ECT (PAT) Sample Cutter (Edge Pressure Sampler)



Equipment Usage

- Products are mainly used for cutting ECT, PAT specimen.

Standard

- GB/T 6546, GB/T6548, GB/T 6544, QB/T 1671 etc.....

Features

- The thrust type sampling, the structure is simple and firm.
- The double stainless steel ball bearing guide sleeve design, equipped with imported Japanese blade demolition convenient, sampling effort, accurate, without burr.

Technical Specifications

Model	JG-349
Maximum sample thickness	18mm
Sample size	100mm, 80mm, 25mm (ECT sample 100 x 25mm; PAT sample 80 x 25)
Sampling accuracy error	±0.5mm
Dimensions	440mm(L) x 380mm (W) x 180mm(H)
Weight	19Kg

