

Thermo Scientific Orion 2100 Series Grab Sampler & FAST QC Module

Rapid analyzer performance verification

Eliminate the guess work and high costs associated with time consuming lab verification of your on-line analyzer - our module for performance verification minimizes calibration frequency with confidence within minutes.



Markets:

- Power
- Pulp and paper
- Chemical/petrochemical

Applications:

- Quality control
- Quality assurance
- Measurement validation and verification
- Grab sample analysis



The new, Thermo Scientific™ Orion™ 2100 series grab sampler and FAST QC™ module provides rapid analysis of analyzer performance with an easy-to-use method. Validate the accuracy of an online analyzer in minutes to ensure optimum performance and response.

Offered as an integrated component of the 2100 series analyzers or as a field upgrade enhancement, the grab sampler and FAST QC module allows maximum flexibility to fit your plant's instrumentation requirements.

Understanding the time and resource constraints for efficient plant operation—the 2100 series calibration and validation process eliminates the need for time consuming laboratory grab sample analysis.

Features

- Offered as integrated module on 2100 series analyzers from the factory or as a field upgradable module at your location
- Rapid delivery of a known standard to verify measurement accuracy—reduces the frequency of time consuming calibrations
- Reaches stability in less than 10 minutes, the system has you back online with confidence

Orion 2100 Series Grab Sampler and FAST QC Module

- Global support – with experience that comes from supporting our online customers for over 35 years throughout the world, our water quality specialists and customer support teams offer a quick, thorough and professional response to any problem.
- Focus on user benefits – we work closely with you to define your needs, and ensure you are using the analyzer in a way that improves your bottom line.

Flexibility for Your Convenience

The 2100 series grab sampler and FAST QC module is offered as integrated component on 2100 series analyzers or as a field upgradable module.

For more information, contact your local water quality specialists, call 1-800-225-1480 or visit www.thermoscientific.com/processwater.



Product Specifications	
Volume	1000 mL
Grab Sample Consumption	500 mL in 10 minutes with an analyzer flow rate of 40 mL/minute
Detection Limits	5 ppb (21GRABLL) 10 ppb (21GRABXP)
Response Time	90% of final reading within 5 minutes after calibration
Compatibility	2100 series analyzers only*
Electrical Requirements	9 VDC, 1 A
Input	115 VAC / 60 Hz adapter, 230 VAC / 50 Hz adapter
Dimensions H x W x D	12.7 cm x 48.3 cm x 34.3 cm (5 inch x 19 inch x 13.5 inch)
Weight	3.6 kg (8 pounds)
Temperature	5 to 45 °C
Approvals	CE, CSA, FCC Class A limits

Ordering Information	
21GRABXP	Grab sampler / FAST QC upgrade module to fit existing 2111XP, 2111XPEN, 2110XP, 2110XPEN, 2117HL, 2117HLEN, 2117XP, 2117XPEN, 2120XP and 2120XPEN analyzers
21GRABLL	Low level sodium grab sampler / FAST QC upgrade module to fit existing 2111LL and 2111LLEN analyzers
21GRBT	Replacement bottle for 2100 series grab sampler
21GRPA	Air pump replacement kit for 2100 series grab sampler
21GRFK	Fittings kit for 2100 series grab sampler

* Not currently available for use with 2109XP, 2109XPEN, 2118XP or 2118XPEN analyzers



Thermo Fisher Scientific
Water Analysis Instruments
Chelmsford, MA USA
Quality Management System
Registered to ISO 9001

thermoscientific.com/processwater

© 2013 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. & its subsidiaries.

Water Analysis Instruments

North America
Toll Free: 1-800-225-1480
Tel: 1-978-232-6000
info.water@thermo.com

Netherlands
Tel: (31) 033-2463887
info.water.uk@thermo.com

India
Tel: (91) 22-4157-8800
wai.asia@thermofisher.com

Japan
Tel: (81) 045-453-9175
wai.asia@thermofisher.com

China
Tel: (86) 21-68654588
wai.asia@thermofisher.com

Singapore
Tel: (65) 6778-6876
wai.asia@thermofisher.com

Australia
Tel: (613) 9757-4300
in Australia (1300) 735-295
InfoWaterAU@thermofisher.com