MYZOX

Multi Cushion Case

Multi Cushion Case M-MCC





Light & Tough

1 Hard Cushion

To protect the instrument from external shock, outer side of the case are surrounded by hard cushions that hardly yield.



2 Free Partition

Removable partition can be attached anywhere by hook and loop tape. Removable softer cushions lead to even higher protection level.

Ideal size for Rotating Laser





3 2Way

Use the handles and the shoulder strap depending on the situation.



4 Multi Case

You can put away various surveying accesories as you like.

Specifications				
Material	Nylon			
Dimensions	Outer : W360 x H250 x D260mm			
	Inner : W320 x H220 x D220mm (without cushions)			
Weight	1.3kg			
Storable model	Myzox : MJ-300, MJ-300S			
	Topcon : RL-H5A, RL-H4C			
	Sokkia : LP510, LP515			
	Nikon : LL300N (1 detector set)			
	Leica : RUGBY610			
	*Please check the inner dimensions for other models.			

Set Contents

■ Soft	cushion	for	the	side	хC	3
■ Soft	cushion	for	the	top	x 1	
■ Shou	Ilder str	an y	1			

- Hard cushion partition x 1 ■ Hard cushion for the bottom x 1
- Handle band x 1

Code No.	Product Name	Model		
222883	Multi Cushion Case	M-MCC		
*All company names and product names used in this document are trademarks				
or registered trademarks of their respective companies.				

Your Local Authorized Myzox Dealer is :

Myzox Co., Ltd.

401 Yamagoe, Nagakute, Aichi 480-1111 Japan TEL: +81-561-64-2210 FAX: +81-561-63-7770

www.myzoxjapan.com