

Oxygen simplified

EverFlo concentrators reduce the cost of oxygen therapy



Ideal oxygen concentrators for you and your patients.

Philips Respironics recognizes the growing need for cost-effective and reliable oxygen concentrators that can reduce the cost of providing oxygen therapy. That's why we've developed a series of forward-thinking oxygen products like EverFlo that simplify your oxygen business.

Expect industry-leading quality and reliability

EverFlo stationary oxygen concentrators can save you time and money over the course of ownership because of their design and quality. Field-proven over more than three years, EverFlo was rated highest in perceived quality by a 2008 independent survey of homecare providers, and hundreds of thousands of units have been deployed worldwide.*

Manufactured in the United States and backed by our deep design knowledge, EverFlo is developed using only high-quality parts. With a recessed flow meter, metal cannula connector, strong plastic casing, proven sieve bed, and Pressure Swing Adsorption (PSA) technology, you will notice a reduction in repairs, complaints, and unscheduled visits.

Enjoy ease of service

EverFlo is designed with fewer parts resulting in lower overall maintenance costs. EverFlo has just two filters. The external filter can be changed in minutes without special skills or tools and the internal micro disk filter is simply checked during routine maintenance procedures.

Differentiate your business

EverFlo concentrators are backed by value-added tools. From in-person and online training to financing options, you have what you need to succeed. And, since Philips Respironics offers a wide range of sleep and home respiratory products and comprehensive support, we can be the source for many of your product needs.

* HME Home Oxygen Survey. 2008. Data on file. Philips Respironics.





Easy to store, transport, and deliver.

EverFlo is compact, has a streamlined and integrated handle, and weighs only 31 pounds. This results in lower shipping and storage costs and may reduce risk of worker injury.



EverFlo is a green product

EverFlo is significantly lighter, quieter, and more energy efficient than its predecessor. These improvements include reductions of weight (40%), energy usage (24%), and packaging (15%).

Green logo

You'll find the Green logo on Philips environmentally friendly products, so that you can make the right environmental choice.



EverFlo is the winner of the Provider's Choice Award at 2008 Medtrade Spring.

Give your patients what they want.

Your patients will appreciate the compact, non-medical design of the EverFlo oxygen concentrator because it will fit right into their home environment. The unit is designed to be small, light, easy to transport and will fit into the trunk of a car. EverFlo has an easy-to-mount humidifier bottle connection and no patient filter.



Molded handle that's easy to pick up and transport with one or both hands



Humidifier bottle holder compatible with all bottle styles

EverFlo product specifications	
Input voltage	120 VAC +/- 10%
Input frequency	60 hz
Average power consumption	350 watts
Oxygen concentration* (at 5 lpm)	93% (+/- 3%)
Liter flow	0.5 to 5 liters per minute
Weight	31 lbs (14 kg)
Sound level	45 dBA typical
Dimensions	23" (584 mm) h, 15" (381 mm) w, 9.5" (241 mm) d
Required maintenance	Compressor intake filter: change every 2 years
OPI (oxygen percentage indicator) alarm levels	Low oxygen: 82%, very low oxygen: 70%
Operating temperature	55°F to 90°F (12°C to 32°C)
Storage/transport humidity	-30°F to 160°F (-34°C to 71°C) up to 95% relative humidity
Operating humidity	Up to 95% relative humidity
Outlet pressure	5.5 PSI
Operating altitude	0 to 7,500' (0 to 2,286 m)

^{*} Device operation above or outside of the specified voltage, LPM, temperature, humidity, and/or altitude values may decrease oxygen concentration levels.

EverFlo ordering information	
Concentrator systems	Part number
EverFlo oxygen concentrator	1020000
EverFlo oxygen concentrator with OPI (oxygen percentage indicator)	1020001
EverFlo Q oxygen concentrator	1020015
EverFlo Q oxygen concentrator with OPI (oxygen percentage indicator)	1020014
Concentrator accessories	Part number
Compressor intake filter	1038831
Micro disk filter	520
Humidifier connector tube	1039642
Stand-alone low-flow flow meter	H537-1
Integrated low-flow flow meter	H644
Flow meter locking kit	528
EverFlo manual in English, French, Spanish, and Portuguese	1020121
EverFlo service manual	1038571
Extended warranty	Part number
Two-year extended warranty with compressor	1038571



Internal tubing storage keeps tubing conveniently available for use



Metal cannula connector for greater durability



Interface panel with on/off switch and alarm indicators



Recessed flow meter reduces the risk of breakage